
O₂-easy

CA-MI Srl

Via Ugo La Malfa 13 - Frazione Pilastro

43013 Langhirano (PR) Italia

Tel. +39 0521 637133 / +39 0521 631138 - FAX +39 0521 639041

E-mail vendite@ca-mi.it / export@ca-mi.it

www.ca-mi.it

Rev. 5 del 19/12/2018

MADE IN CHINA

O₂-easy

- IT Manuale d'uso
- EN Instruction manual
- DE Handbuch
- FR Mode d'emploi
- ES Manual de instrucciones
- PT Manual de uso
- PL Instrukcja obsługi
- RU Инструкция по эксплуатации

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Thank you for choosing our O2-EASY device for measuring oxygen saturation in the blood.

Oxygen saturation is a percentage of the amount of oxyhaemoglobin (HbO₂) - oxygen combined with an amount of haemoglobin (Hb) - in the blood. In other words, it represents the consistency of oxyhaemoglobin in the blood, a very important parameter for the respiratory system. Many respiratory diseases may cause lowering of the oxygen saturation threshold in human blood. Moreover, the following factors may reduce oxygen saturation: automatic regulation of the dysfunction of an organ caused by anaesthesia, intensive post-surgery trauma, and injury caused by medical examinations. All this may cause a feeling of weakness, asthenia and nausea. It is therefore very important to know what the oxygen saturation of a patient is, so that the doctor can promptly solve any related problems.

The small finger pulse oximeter consumes little power and is easy to use and carry around.

For diagnosis, all the patient has to do is insert the tip of a finger in the photoelectric sensor and the oxygen saturation value will be shown on the display.

Clinical trials have shown the high precision and repeatability of this device.

PARTS DESCRIPTION

- 1 Cord
- 2 Two AAA batteries
- 3 O2-Easy

MEASURING PRINCIPLES

The principle at the heart of the pulse oximeter is a mathematical formula, established using the Lambert Beer Law according to the characteristics of the Absorption Spectrum of Reduced Haemoglobin (RHb) and Oxyhaemoglobin (HbO₂) in light and infrared areas.

The technology adopted is photoelectric inspection and recording of the oxyhaemoglobin according to the pulse capacity, so that two light beams of different wavelength (660nm and 905nm) can be focused on a human fingertip through a blocking sensor. The measurement obtained from the light-sensitive element is shown on the display by means of a process in the electronic circuits and a microprocessor.

OPERATING DIAGRAM

- 1 Infrared ray emission tube
- 2 Infrared ray reception tube

CAREFULLY READ THE USER MANUAL BEFORE EACH USE.

NEVER DISMANTLE THE DEVICE. FOR ANY OPERATION ON THE DEVICE,
CONTACT CA-MII TECHNICAL SERVICE.

THE DEVICE IS NOT SUITABLE FOR CONTINUOUS MONITORING.

SUITABLE FOR THE DETECTION ON CHILDRENS
OF AGE NOT LOWER THAN 3 YEARS.

PRODUCT FEATURES

- 1 The device is easy and practical to use. The device is small, lightweight and easy to carry around.
- 2 The device has low energy consumption and the two AAA batteries have a life of 30 hours of continuous use.
- 3 When the battery charge is so low that it may affect normal functioning of the pulse oximeter, a warning is shown on the display.
- 4 The device automatically turns off when it does not receive any signals for more than 8 seconds.

FIELD OF APPLICATION

The Fingertip Pulse Oximeter is a handheld non-invasive device intended for spot-checking of oxygen saturation of arterial hemoglobin (SpO₂) and Pulse Rate of adult, adolescent and child patients at home (HOME-CARE) or in hospital (including clinical use in surgery, anaesthesia, intensive care, etc.).

The device is not intended for continuous monitoring.

The pulse oximeter does not require any routine calibration or other maintenance operations apart from replacing the batteries. The device is not intended for use during patient transport outside the health facility.

GENERAL SAFETY RULES

- 1 Operation of the fingertip pulse oximeter may be affected by the use of an electrosurgical unit (ESU).
- 2 Do not use the pulse oximeter when the patient is undergoing MRI and CT analysis.
- 3 Do not use the fingertip pulse oximeter in situations where alarms are required. The device has no alarms. It is not for continuous monitoring.
- 4 The fingertip pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
- 5 Frequently check that the pulse oximeter sensor is positioned correctly in relation to the circulation and sensitivity of the patient's skin
- 6 Prolonged use or the patient's conditions might require periodically changing the position of the sensor. Change the position of the sensor and check the integrity of the skin, the circulation and proper alignment at least every 4 hours.
- 7 Autoclave sterilization processes, ethylene oxide sterilization, or immersing the sensor in liquid substances may cause inaccurate measurements or damages.
- 8 Never immerse the device in water.
- 9 Do not use the device in places where there are anesthetic mixtures flammable in contact with air, oxygen or nitrous oxide.

- 10 Do not touch the device with wet hands and never let the device come into contact with liquids.
- 11 If the device is used by children and/or incapable persons, they must always be closely supervised by an adult with full mental faculties.
- 12 Store and use the device in environments protected against atmospheric agents and away from any heat sources.
- 13 Any repair operations must exclusively be carried out by Manufacturer/Distributor Technical Service or by a technical service authorized by the manufacturer using original spare parts. Failure to comply with the above listed rules may compromise the safety of the device.
- 14 This device must be used exclusively for the purpose for which it was designed as described in this manual.
- 15 If the device is not used for an extended period of time, it is advisable to remove the batteries.
- 16 Using the device in ambient conditions different from those indicated in this manual may compromise its safety and parameters.
- 17 This device can also be used on children older than 3 years, provided that there is a normal size-weight / age ratio of the limb where it's applied.

WARNING!

Below are the possible conditions / uses that can cause a misreading or inaccurate measurement of the device:

- 1 Significant levels of dysfunctional hemoglobin (such as carbonyl - hemoglobin or methemoglobin).
- 2 Intravascular dyes such as indocyanine green or methylene blue.
- 3 SpO₂ measurements may be negatively affected where there is too much ambient light. Cover the sensor if necessary (for example, with a surgical towel).
- 4 Excessive patient movement.
- 5 High-frequency electrosurgical interference and defibrillators.
- 6 Venous pulsations.
- 7 Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
- 8 The patient has hypertension, severe vasoconstriction, severe anemia, or hypothermia.
- 9 The patient is in cardiac arrest or is in shock.

EN

- 10 Fingernail polish or false fingernails.
- 11 Weak pulse quality (low perfusion).
- 12 Low hemoglobin.

Carefully read these instructions before using the device. Not following these instructions may cause incorrect measurements or damage to the device or the user.

The manufacturer is not responsible for the safety, reliability and proper functioning of the device nor for monitoring errors or damage to persons and the product if the user fails to read these instructions, and such events are not covered by the warranty.

TECHNICAL SPECIFICATIONS

Display type:	LED
SpO ₂ :	Measuring interval: 0-100% Accuracy: 70%-100%, ±2%; ≤69% no definition Resolution: 1%

31

EN

Pulse frequency:	Measuring interval: 30-250 BPM Accuracy: 30~99bpm, ±2bpm; 100~250bpm, ±2% Resolution: 1bpm Pulse intensity: Bar graph indicator
Power supply requirements:	Two AAA alkaline batteries 1.5V - Lifetime until 30 hours Energy consumption: Less than 40mA Battery charge indication: 
Dimensions:	Length: 57mm~59mm Width: 33mm~35mm Height: 34mm~36mm
Weight:	40g~60g (including two AAA batteries)
Ambient requirements:	Operating temperature: 5~40°C Storage temperature: -25~70°C Ambient humidity: 15~93%, without formation of condensate during operation. ≤93%, without formation of condensate during storage. Atmospheric pressure: 70kPa~106kPa

32

Equipment data update period:
Data update period of slower average is 8s (Fig A).


Fig A

ARMS Value Analysis Statement

Item	90-100	80-90	70-80
#pts	78	66	63
Bias	1.02	0.40	-0.48
ARMS	1.66	1.46	1.63


Fig B

33

The following details are provided to disclose actual performance observed in the clinical validation study of healthy adult volunteers. The ARMS value analysis statement and Bland-Altman plot of data is shown as in the Fig B picture.


MEASURING PROCEDURE

- 1 Fit two AAA batteries in the battery compartment and close the cover.
- 2 Open the clamp as shown in the figure.
- 3 Insert a finger in the pulse oximeter.
- 4 Press the button on the front panel once.
- 5 The measurement is shown on the display.


Note:

Only Use alcohol to clean the rubber part of pulse oximeter in contact with the patient's finger (the rubber part in the device is medical rubber that does not contain toxins and is not harmful to the skin).

Insert the finger in the pulse oximeter with the nail facing up.


34

There are two display modes:

After turning on the pulse oximeter, each time you press the power switch, the pulse oximeter will switch to another display modes.

The bar graph displays the pulse measured on the patient's wrist. The height of the bar graph shows the pulse resistance in the patient's wrist.


DEVICE ACCESSORIES

- 1 Cord
- 2 Two AAA batteries
- 3 Instruction manual


CORD INSTALLATION

- A Slip the thin part of the cord through the eyelet.
- B Position the end into place before pulling.

BATTERY INSTALLATION

- 1 Fit the two AAA batteries provided in the battery compartment respecting the polarity.
- 2 Close the cover of the battery compartment horizontally along the arrow.

**Notes:**

- The batteries must be fitted respecting the correct polarity, otherwise the device may be damaged.
- Fit and remove the batteries in the right order to prevent damaging the device.
- Remove the batteries if the pulse oximeter is not used for a long period of time.

MAINTENANCE AND STORAGE CONDITIONS

- 1 Replace the batteries when the warning light comes on.
- 2 Clean the surface of the pulse oximeter before using it for diagnosis.
- 3 Remove the batteries from the battery compartment if the device is not used for a long period of time.
- 4 It is advisable to store the device at a temperature between -25°C and 70°C with relative humidity $\leq 93\%$.
- 5 The device must always be stored in a dry environment. A humid environment may damage the device reducing its life.
- 6 Follow the Italian regulations for product and spent battery disposal.

CLEANING THE FINGERTIP PULSE OXIMETER

Please use medical alcohol to clean the silicone touching the finger inside of oximeter with a soft cloth dampened with 70% isopropyl alcohol. Also clean the being tested finger using alcohol before and after each test. Do not pour or spray liquids onto the oximeter, and do not allow any liquid to enter any openings in the device. Allow the oximeter to dry thoroughly before reuse.

The fingertip pulse oximeter requires no routine calibration or maintenance other than replacement of batteries.

The use life of the device is five years when it is used for 15 measurements every day and 10 minutes per one measurement.

Stop using and contact local service center if one of the following cases occurs:

- An error in the Possible Problems and solutions is displayed on screen.
- The oximeter cannot be powered on in any case and not the reasons of battery.
- There is a crack on the oximeter or damage on the display resulting readings cannot be identified; the spring is invalid; or the key is unresponsive or unavailable.

DISINFECTING:

The applied parts touching the patients' body are required to be disinfected once after each use. The recommended disinfectants include: ethanol 70%, isopropanol 70%, glutaraldehyde-type 2% liquid disinfectants.

Disinfection may cause damage to the equipment and is therefore not recommended for this pulse oximeter unless otherwise indicated in your hospital's servicing schedule. Clean the pulse oximeter before disinfecting it.

CAUTION: Never use EtO or formaldehyde for disinfection.

TROUBLESHOOTING

Problem	Possible causes	Solutions
The SpO2 or PR values are not normally displayed	<ol style="list-style-type: none"> 1 The finger is incorrectly inserted. 2 The patient's oxyhaemoglobin value is too low to be measured. 	<ol style="list-style-type: none"> 1 Insert the finger correctly and retry. 2 Try a few more times if you are sure that the problem is not due to device malfunctioning. See a doctor for an exact diagnosis.
The SpO2 or PR value are not stably shown	<ol style="list-style-type: none"> 1 The finger might not be inserted deep enough 2 The finger is unstable or the patient has moved. 	<ol style="list-style-type: none"> 1 Insert the finger correctly and retry. 2 Retry without moving.

TROUBLESHOOTING

Problem	Possible causes	Solutions
The pulse oximeter does not turn on	<ol style="list-style-type: none"> 1 The batteries might be flat or not fitted. 2 The batteries might be incorrectly fitted. 3 The pulse oximeter might be damaged. 	<ol style="list-style-type: none"> 1 Replace the batteries. 2 Remove and refit the batteries. 3 Contact your local Customer Service.
The warning lights unexpectedly go off	<ol style="list-style-type: none"> 1 The device is programmed to automatically turn off when it does not receive a signal for more than 8 seconds. 2 The batteries are low. 	<ol style="list-style-type: none"> 1 No problem, this condition is normal. 2 Replace the batteries.
"Error7" appears on the display	<ol style="list-style-type: none"> 1 Batteries low. 2 The emission tube is damaged. 3 The control circuit is malfunctioning. 	<ol style="list-style-type: none"> 1 Replace the batteries. 2 Contact your local Customer Service. 3 Contact your local Customer Service.

SYMBOLS

	Type BF applied part
	Attention
%SpO ₂	Oxygen saturation
	Pulse rate (BPM)
	Low power indication
SN	Serial No.
	No SpO ₂ alarm
	Power switch
CE 0123	Mark of compliance to the MDD 93/42 EEC
	Follow instruction for use
	Manufacturer: CA-MI Srl – Via Ugo La Malfa 13 Frazione Pilastro – 43013 Langhirano (PR) Italy
	Storage relative humidity
	Storage temperature

SYMBOLS

IP22	Protection rating under solid or liquid penetration of electronic equipment or to an enclosure for electronic equipment.	
	1° NUMBER Solid penetration. Protected against solid objects over Ø 12mm e.g. hands, large tools.	2° NUMBER Liquid penetration. Protected against falling drops of water, if the case is disposed up to 15° from vertical.

- This equipment complies with IEC 60601-1-2:2014 for electromagnetic compatibility for medical electrical equipment and/or systems. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare and other environments, it is possible that high levels of such interference due to close proximity or strength of a source might disrupt the performance of this device.

- These product's materials that contact with the patient's skin contain medical silicone and ABS plastic enclosure are all pass the ISO10993-5 Tests for invitro cytotoxicity and ISO10993-10 Tests for irritation and delayed-type hypersensitivity.

IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2012/19/UE-WEEE:

The symbol on the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, don't dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the sample type to be used with the same functions. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same. Caution: The wrong disposal of electric and electronic equipment may involve sanctions.

DISPOSAL OF WASTE BATTERIES (Directive 2006/66/EC):

This symbol on the battery or on the packaging indicates that the battery provided with this product shall not be treated as household waste. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. The recycling of the materials will help to conserve natural resources. At the end of their life hand the batteries over to the applicable collection points for the recycling of waste batteries. For more detailed information about recycling of this product or battery, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

EN

RISK OF ELECTROMAGNETIC INTERFERENCE AND POSSIBLE REMEDIES:

The O2-EASY is an electro-medical device that requires particular precautions regarding electro-magnetic compatibility and which must be installed and commissioned according to the electro-magnetic compatibility information supplied.

Portable and mobile radio communication devices (mobile phones, transceivers, etc.) may interfere with the medical device and should not be used in close proximity with, adjacent to or on top of the medical device. If such use is necessary and unavoidable, special precautions should be taken so that the electro-medical device functions properly in its intended operating configuration (for example, constantly and visually checking for the absence of anomalies or malfunctions). The use of accessories, transducers and cables different to those specified, with the exception of transducers and cables sold by the appliance and system manufacturer as spare parts, can lead to an increase in emissions or in a decrease of the immunity of the device or system.

45

FACTORY WARRANTY

CA-MI Srl warrants its products for 24 months after purchasing date. In front of this warranty, CA-MI Srl will be obliged only to repair or substitute free of charge the products or parts of them that, after verification effected on our factory, or our authorized Service Center, by the Technical Service, results defective.

CA-MI Srl will substitute the device, only if a SALE RECEIPT. We also request to specify the kind of fault, in order to speed up repairing procedures. The product must be accompanied by a description of the defect. The warranty, with exclusion of responsibility for direct and indirect damages, is thought limited to the solos defects of material or workmanship and it stops having effect when the device results however gotten off, tampered or sheltered out of the Factory or from the Authorized Service center.

The commodity always travels to risk and danger of the buyer, without any responsibility of CA-MI Srl for damages caused by the transport or dismay from the vector.

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