

MATTRESS



premium hospital beds



1946-2016
70 years

IMO



medical beds worldwide



CORPORATE PROFILE

- COMPANY **IMO-Indústrias Metalúrgicas, S.A.**
- HEAD-QUARTERS **Longra, Felgueiras, Portugal**
- VAT No. **500.136.939**
- FOUNDATION YEAR **1946**
- EMPLOYEES **112**
- SHARE CAPITAL **2,5 M€**
- TURN-OVER (2015) **5,7 M€**
- INSTALLED PRODUCTS **46 países**
- INDUSTRIAL FACILITIES **8.400 m2 área coberta**
- PRODUCTION CAPABILITIES
 - Research and development**
 - Prototyping and tools**
 - Molds Engineering**
 - Laboratory quality control**
 - Electronics**
 - Metalworking**
 - Plastics**
 - Woodwork**
 - Upholstery**
 - 2 assembly lines**
 - Logistics and distribution**


FLUXAIR



DENSITY	30 Kg/m ³ ±2,0
RESISTANCE TO TRACTION	Min. 1,30 Kpa
RESILIENCE	Min. 35 %
PERMANENT DEFORMATION	8% (A 75% Max.)
FATIGUE TEST (Hardness)	Max. 35%
FATIGUE TEST (Height)	Max. 2.0 %

fluxair mattress

Mattress made in a polyurethane foam body sectioned as cubes mesh. This mesh creates inner cavities that contains the necessary air to the mattress cover unidirectional ventilation effect. Those cubes distribution does a pressure relief alternation to get an anti-decubitus effect. It includes a Telastic sanitary cover.



RUPTURE FORCE	>600 N/5cm
RUPTURE STRETCH	>60%
WATER-PROOF	OK (ISO811-1981)
LAUNDRY SCHRINK	0% (ISO5077-1984)
FLAMMABILITY	OK Crib5
BREATHABILITY 37°/65%r.h.	>1000 g/m (24h)

telastic sanitary cover

Made in a very thin, soft and elastic polyurethane fabric, the cover copy perfectly the bed articulations. Is fire retardant and water proof . The fabric mesh is made by a conic shape performing a breathable effect. Its composition has a antialergic treatment that prevents the microorganisms spread.

IMO-INDUSTRIAS METALURGICAS, S.A.

www.imo.com.pt

HEADQUARTERS & FACTORY

Rua Dr. Aurélio Teixeira de Sousa, 575
4650-312 Rande
PORTUGAL
T. +351 255 340 220
F. +351 255 340 239
sede@imo.com.pt

IMO PORTO

Rua Adolfo Casais Monteiro, 117
4050-520 Porto
Tel.: 222 009 767
Fax: 222 000 647
porto@imo.com.pt

IMO LISBOA

Av. Defensores de Chaves, 87-A
1000-115 Lisboa
Tel.: 217 963 478
Fax.: 217 975 822
lisboa@imo.com.pt

AUTHORIZED DISTRIBUTOR: