

## Q500® M/F/R with SEDEO® PRO



## ORDER FORM


Errors and imissions excepted.

Billing Address Or Customer-ID			Shipping Address	Identical To Billing Address
C.-ID:			C.-ID:	
Name:			Name:	
Street:			Street:	
Town/PC:			Town/PC:	
E-Mail:			E-Mail:	
Phone:			Phone:	
Order Information			User Weight	
Order	Calculation	Built-4-Me®	<b>User Weight Capacity: 160 kg</b>	
Purchase Order No.:			<b>User Weight:</b> kg	
<b>Art. No.:</b> Q500 M: 2250050M-000, Q500 F: 2250050F-000, Q500 R: 2250050R-000			Reimbursement Code: xxx	
Order Date:	Day	Month	Year	

Base Price € **XX.XXX,-**

## Client Measurements

a) Top Of Head To Seat	mm
b) Top Of Shoulder To Seat	mm
c) Underarm To Seat	mm
d) Elbow To Seat	mm
e) Top Of Pelvis To Seat	mm
f) Shoulder Width	mm
g) Trunk Width	mm
i) Hip Width	mm
j) Actual Sitting Depth	mm
m) Foot Length	mm
n) Lower Leg Length	mm
o) Actual Seat Width	mm



## Seating Specifications

All Bases Prepared For DAHL Docking Station		<b>Seat Type</b>	SEDEO PRO
<b>Typical WC Weight</b>	Without Batteries 140 kg	<b>Seat Size Range</b>	Width 380 - 560 mm, Depth 400 - 560 mm
	With Batteries 165 kg	<b>Seat To Floor Height Range</b>	415 mm - 490 mm

## MWD Base Specifications

<b>User Weight Capacity</b>	160 kg	<b>Drive Wheels</b>	13" Drive Wheels (Black)
<b>Base Size</b>	Width Min. 620 mm	<b>Tires / Castors</b>	7" Front / 7" Rear (Black)
<b>Drive Control</b>	VR2, Option R-Net	<b>Charger</b>	8 A Dual Mode Charger
<b>Driving Range</b>	Up To 40 km	<b>Battery Size</b>	60 And 80 Ah Batteries
<b>Drive Motors</b>	4-Pole 6 kph, Options 8, 10 And 12.5 kph	<b>Colour Accents</b>	7 Colour Accents
<b>Ground Clearance</b>	70 mm	<b>Frame Colour</b>	Black

## FWD Base Specifications

<b>User Weight Capacity</b>	160 kg	<b>Drive Wheels</b>	13" / 14" Drive Wheels (Black)
<b>Base Size</b>	Width Min. 610 mm	<b>Tires / Castors</b>	9" Rear (Black), Options 10"
<b>Drive Control</b>	VR2, Option R-net	<b>Charger</b>	8 A Dual Mode Charger
<b>Driving Range</b>	Up To 40 km	<b>Battery Size</b>	60 And 80 Ah Batteries
<b>Drive Motors</b>	4-Pole 6 kph, Options 8, 10 And 12.5 kph	<b>Colour Accents</b>	7 Colour Accents
<b>Ground Clearance</b>	65 mm	<b>Frame Colour</b>	Black

## RWD Base Specifications

<b>User Weight Capacity</b>	160 kg	<b>Drive Wheels</b>	13" / 14" Drive Wheels (Black)
<b>Base Size</b>	Width Min. 610 mm	<b>Tires / Castors</b>	9" Front (Black), Options 10"
<b>Drive Control</b>	VR2, Option R-net	<b>Charger</b>	8 A Dual Mode Charger
<b>Driving Range</b>	Up To 40 km	<b>Battery Size</b>	60 And 80 Ah Batteries
<b>Drive Motors</b>	4-Pole 6 kph, Options 8, 10 And 12.5 kph	<b>Colour Accents</b>	7 Colour Accents
<b>Ground Clearance</b>	65 mm	<b>Frame Colour</b>	Black

## Base

## Q500 Drive Bases

€

QOP000023	<b>Q500 M (MWD Base Chair)- Should be order through UK*</b>	XXX,-
QOP000024	<b>Q500 F (FWD Base Chair)</b>	XXX,-
QOP000029	<b>Q500 R (RWD Base Chair)</b>	XXX,-

## Motor Packages

€

## MWD Base

QOP180011	2-Pole Motor, 6 kph/ 4 mph, Standard	XXX,-
-----------	--------------------------------------	-------

## MWD/FWD/RWD Base

<input type="checkbox"/>	QOP180012	2-Pole Motor, 6kph/ 4mph, Performance	XXX,-
	QOP180009	4-Pole Motor, 6 kph / 4 mph (8 kph Motor With Reduced Speed)	XXX,-
	QOP180010	4-Pole Motor, 8 kph / 5 mph	XXX,-
	QOP180007	4-Pole Motor, 10 kph / 6.5 mph (1) (2)	XXX,-
	QOP180008	4-Pole Motor, 12.5 kph / 8 mph (1) (3)	XXX,-

## Tracking

€

QOP110009	Traction Control System For R-Net (GYRO Technology)	XXX,-
QOP110014	Omit Traction Control System For R-Net	XXX,-

Drive Wheels 

€

## MWD/FWD/RWD Base

QOP070065	Drive Wheel, 13", Pneumatic, Knobbly, Black	XXX,-
QOP070069	Drive Wheel, 13", Solid, Knobbly, Black	XXX,-

## FWD/RWD Base

QOP070032	Drive Wheel, 14", Pneumatic, V-Type, Black	XXX,-
QOP070031	Drive Wheel, 14", Solid, V-Type, Black	XXX,-
QOP070034	Drive Wheel, 14", Pneumatic, Knobbly, Black	XXX,-
QOP070033	Drive Wheel, 14", Solid, Knobbly, Black	XXX,-

Castor Wheels 

€

## MWD Base

QOP080213	Castor Wheel, 7", Solid, Black	XXX,-
-----------	--------------------------------	-------

## FWD Base

QOP080210	Castor Wheel, 9", Solid, Black	XXX,-
QOP080211	Castor Wheel, 9", Pneumatic, Black	XXX,-

## RWD Base

QOP080210	Castor Wheel, 9", Solid, Black	XXX,-
QOP080211	Castor Wheel, 9", Pneumatic, Black	XXX,-
QOP080212	Castor Wheel, 10", Pneumatic, Black	XXX,-

Mudguards (FWD/RWD Only) 

€

QOP010007	Mudguard, Drive Wheel, 14"	XXX,-
-----------	----------------------------	-------

(1) Highly Recommended With 80 Ah Battery

(2) If FWD Selected: Highly Recommended With Traction Control System

(3) If FWD Selected: Only In Combination With Traction Control System

QOP010008	Mudguard, Castor Wheel, 9"	XXX,-
-----------	----------------------------	-------

### Kerb Climber (RVD Only)

€

QOP090038	Kerb Climber	XXX,-
-----------	--------------	-------

### Base

### Colour Accents

€

QOP140005	Apple Green (RAL 110 70 77)	XXX,-
QOP140004	Black (RAL 9005)	XXX,-
QOP140003	Blue (RAL 260 40 45)	XXX,-
QOP140006	Orange (RAL 060 60 80)	XXX,-
QOP140002	Red (RAL 040 40 67)	XXX,-
QOP140007	Silver (RAL 9006)	XXX,-
QOP140001	White (RAL 9010)	XXX,-

### LED Lights & Indicators

€

QOP110201	LED Light And Indicator StVZO Approved (1)	XXX,-
QOP110209	Omit Light And Indicator	XXX,-

### Seat Heights (2)

€

SH415	415 mm (3) (4)	-
SH425	425 mm (5)	
SH44	440 mm (4)	-
SH45	450 mm (5)	
SH465	465 mm	-
SH49	490 mm	-

### Batteries & Chargers

€

#### Batteries

QOP160003	Battery, 50Ah, Gel	XXX,-
QOP160001	Battery, 60 Ah, Gel, Maintenance-Free	XXX,-
QOP160005	Battery, 80 Ah, Gel, Maintenance-Free (6)	XXX,-
QOP160009	Omit Batteries	XXX,-

#### Chargers

QOP160101	Charger, Automatic, 220 V / 8 A	XXX,-
QOP160100	Charger, Automatic, 220 V / 8 A, UK	XXX,-
QOP160102	Charger, Classic, 220 V / 8 A, IP54	XXX,-
QOP160103	Charger, Automatic, 220 V / 10 A, UK	XXX,-
QOP160104	Charger, Automatic, 220 V / 10 A	XXX,-
QOP160109	Omit Charger	XXX,-

(1) QOP110030 Lighting Module VR2 Or QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required

(2) Measured From Front Of Seat Pan To Floor. Obtain Quote To Verify Seat To Floor Heights

(3) Only With 60 Ah Batteries

(4) Not With QOP130009 Seat Elevator/Tilt Module, Powered

(5) Only For QOP130009 Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered

(6) Highly Recommended With Battery Charger 220 V / 10 A

## Powered Seating

### Power Seating Options

€

QOP130011	Fixed Tilt Module (Adjustable 0°, 3°, 6°, 9°, Factory Setting 6°)	–
QOP130004	Posterior Tilt, Powered, 0° - 30° (1)	–
QOP130012	Posterior Tilt, Powered, -5° - 25° (1) (2)	–
QOP130007	Posterior Tilt, Powered, 0° - 50° (1)	–
QOP130005	Posterior Tilt, Powered, -5° - 45° (1) (2)	–
QOP130009	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered (1)	–

### Backrest Frame Type

€

QOP030083	Backrest, Tool Angle Adjustable, 90° - 120°	–
QOP030084	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Right (3)	XXX,-
QOP030085	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Left (3)	XXX,-
QOP030086	Backrest, Powered Recline, 90° - 130° (1)	XXX,-
QOP030087	Backrest, Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear (1)	XXX,-

### Backrest Frame Option

€

QOP030098	Backrest, Flip Down, Quick Release	XXX,-
-----------	------------------------------------	-------

### Back Cushion Type

€

#### SEDEO PRO (See Section A)

QOP030104	SEDEO PRO, Back Cushion	–
-----------	-------------------------	---

#### Interfaces

QOP030106	SEDEO PRO, Interface Only	–
QOP030107	SEDEO PRO, Interface Only For Back Cushion, Flexible	–
QOP030209	Omit Back Cushion	–

### A – SEDEO PRO Back Size

€

#### Back Cushion Width

BW38	380 mm, For Effective Back Width 380 mm	–
BW42	420 mm, For Effective Back Width 380 - 420 mm	–
BW46	460 mm, For Effective Back Width 420 - 460 mm	–
BW50	500 mm, For Effective Back Width 460 - 500 mm	–
BW54	540 mm, For Effective Back Width 500 - 540 mm	–
BW56	560 mm, For Effective Back Width 520 - 560 mm	–

#### Back Cushion Height

BH44	440 mm, For Effective Back Height 440 - 500 mm (4)	–
BH50	500 mm, For Effective Back Height 500 - 560 mm	–
BH56	560 mm, For Effective Back Height 560 - 620 mm	–

### A – SEDEO PRO Back Cushion

€

QOP030109	SEDEO PRO, Back Cushion, Profiled	XXX,-
QOP030113	SEDEO PRO, Back Cushion, Profiled, Adjustable Straps	XXX,-
QOP030114	SEDEO PRO, Back Cushion, Deep Contour, Passive	XXX,-
QOP030115	SEDEO PRO, Back Cushion, Flexible, Deep Contour, Passive, Adjustable Straps	XXX,-
QOP030116	SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps	XXX,-
QOP030117	SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps	XXX,-
	SEDEO PRO, Back Cushion, Adaptive (Choose B4Me Requirements Separately)	XXX,-

(1) Not With QOP110036 Remote Joystick, VR2. In Case Of R-Net: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L

(2) Not With Seat Height 415 Or 440 mm

(3) Not On Same Joystick Side

(4) Only With QOP030116 SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps & QOP030117 SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps

## Powered Seating

### A – SEDEO PRO Back Cushion Cover

€

QOP030059	Back Cushion Cover, Incontinence Protection	–
QOP030108	Back Cushion Cover, Comfort	XXX,-

### A – SEDEO Lateral Thoracic Supports

€

QOP030607	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Right	XXX,-
QOP030608	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Left	XXX,-
QOP030609	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Short Pad, 90 x 120 mm, Right	XXX,-
QOP030610	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Short Pad, 90 x 120 mm, Left	XXX,-

### Headrest Pad

€

QOP030430	SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm (1)	XXX,-
QOP030431	SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm (1)	XXX,-
QOP030432	SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm (1)	XXX,-
WHT030401	Whitmyer, Headrest Pad, Plush, (LxH) 203 x 127 mm (2)	XXX,-
WHT030402	Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm (2)	XXX,-
WHT030403	Whitmyer, Headrest Pad, Plush, (LxH) 356 x 127 mm (2)	XXX,-
WHT030404	Whitmyer, Headrest Pad, Adjust-A-Plush, Small, (LxH) 355 x 76 mm (2)	XXX,-
WHT030406	Whitmyer, Headrest Pad, Adjust-A-Plush, Large, (LxH) 460 mm x 127 mm (2)	XXX,-

### Headrest Mount (3)

€

QOP030433	SEDEO, Headrest Mount	XXX,-
QOP030434	SEDEO, Headrest Mount, With Quick Release	XXX,-
WHT030407	Whitmyer, Headrest Mount, LINX	XXX,-
WHT030408	Whitmyer, Headrest Mount, LINX <sup>2</sup>	XXX,-
WHT030409	Whitmyer, Headrest Mount, AXYS	XXX,-
WHT030410	Whitmyer, Headrest Mount, AXYS, Offset	XXX,-
WHT030411	Whitmyer, Headrest Mount, COBRA XTRA Flip Back	XXX,-

### Headrest Receiver

€

QOP030119	Headrest Receiver, Back Cushion Mount (4)	XXX,-
QOP030120	Headrest Receiver, Backrest Frame Mount	XXX,-

### Push Handles

€

QOP030204	Push Handles, Height-Adjustable	XXX,-
-----------	---------------------------------	-------

### Seat Rail

€

QOP020085	Seat Rail, Short, For Seat Depth 400 - 520 mm	–
QOP020069	Seat Rail, Long, For Seat Depth 460 - 560 mm (5)	XXX,-

(1) Only Use With SEDEO Headrest Mount

(2) Only Use With Whitmyer Headrest Mount

(3) Choose A Whitmyer Headrest Mount And Headrest Receiver In Case You Would Like To Add A Head Control In The Special Control Section

(4) Only With QOP030109 SEDEO PRO, Back Cushion, Profiled &amp; QOP030113 SEDEO PRO, Back Cushion, Profiled, Adjustable Straps &amp; QOP030114 SEDEO PRO, Back Cushion, Deep Contour, Passive

(5) Not With Seat Width 380 &amp; 400 mm

## Powered Seating

### Seat Cushion Type

€

#### SEDEO (See Section C)

QOP020027	SEDEO, Seat Cushion	-
-----------	---------------------	---

#### JAY (See Section D)

QOP020515	JAY, Seat Cushion	XXX,-
-----------	-------------------	-------

#### Omit (Select Seat Width And Depth in Section C)

QOP020070	Omit Seat Cushion	XXX,-
-----------	-------------------	-------

### C – SEDEO Seat Size

€

#### Seat Cushion Width

SW38	380 mm, For Effective Seat Width 380 mm	-
SW42	420 mm, For Effective Seat Width 380 - 420 mm	-
SW46	460 mm, For Effective Seat Width 420 - 460 mm	-
SW50	500 mm, For Effective Seat Width 460 - 500 mm	-
SW54	540 mm, For Effective Seat Width 500 - 540 mm	-
SW56	560 mm, For Effective Seat Width 520 - 560 mm	-

#### Seat Cushion Depth

SD40	400 mm, Comes With Cushion Size 460 mm	-
SD42	420 mm, Comes With Cushion Size 460 mm	-
SD44	440 mm, Comes With Cushion Size 460 mm	-
SD46	460 mm, Comes With Cushion Size 500 mm	-
SD48	480 mm, Comes With Cushion Size 500 mm	-
SD50	500 mm, Comes With Cushion Size 500 mm	-
SD52	520 mm, Comes With Cushion Size 560 mm	-
SD54	540 mm, Comes With Cushion Size 560 mm	-
SD56	560 mm, Comes With Cushion Size 560 mm	-

### C – SEDEO Seat Cushions

€

QOP020073	SEDEO, Seat Cushion, Profiled, Height 70 mm	-
QOP020074	SEDEO, Seat Cushion, Memory, Height 70 mm (1)	XXX
QOP020075	SEDEO, Seat Cushion, Support, Height 90 mm	XXX
QOP020076	SEDEO, Seat Cushion, Support Memory, Height 90 mm	XXX
QOP020077	SEDEO, Seat Cushion, Form, Height 110 mm (1)	XXX
QOP020078	SEDEO, Seat Cushion, Stabilizer, Height 90 mm	XXX
	SEDEO, Seat Cushion, Adaptive (Choose B4Me Requirements Separately)	XXX

### C – SEDEO Seat Cushion Cover

€

QOP020518	Seat Cushion Cover, Incontinence Protection	XXX
QOP020080	Seat Cushion Cover, Comfort	XXX

### C – SEDEO Seat Boards

€

QOP020082	SEDEO, Seat Board, Flat, Depth-Adjustable, Cushions Guide At Rear	XXX
QOP020083	SEDEO, Seat Pan, 50 mm, Depth-Adjustable, Cushions Guide At Rear	XXX
QOP020084	SEDEO, Seat Pan, 75 mm, Depth-Adjustable, Cushions Guide At Rear (2)	XXX

(1) Only With QOP020518 Seat Cushion Cover, Incontinence Protection, Only For NL

(2) Only For NL

## Powered Seating

### D – JAY Seat Size

€

#### Seat Cushion Width

SW40	400 mm	–
SW42	420 mm	–
SW44	440 mm	–
SW46	460 mm	–
SW48	480 mm	–
SW50	500 mm, Comes With Cushion Size 500 mm	–
SW52	520 mm, Comes With Cushion Size 500 mm	–
SW54	540 mm, Comes With Cushion Size 560 mm	–
SW56	560 mm, Comes With Cushion Size 560 mm	–

#### Seat Cushion Depth

SD40	400 mm	–
SD42	420 mm	–
SD44	440 mm	–
SD46	460 mm	–
SD48	480 mm	–
SD50	500 mm	–
SD52	520 mm	–
SD54	540 mm (1)	–
SD56	560 mm (1)	–

### D – JAY Seat Cushions

€

JAY020005	JAY Easy Fluid With Flat Base	XXX,-
JAY020003	JAY Soft Combi P	XXX,-
JAY020006	JAY Easy Visco With Flat Base	XXX,-
JAY020010	JAY Balance, Air, Single Valve, Microclimatic Cover, Until Seat Width 500 mm	XXX,-
JAY020012	JAY Balance, Air, Single Valve, Microclimatic Cover, From Seat Width 560 mm	XXX,-
JAY020014	JAY Balance, Air, Single Valve, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm	XXX,-
JAY020016	JAY Balance, Air, Single Valve, Microclimatic Cover, With Positioning Kit, From Seat Width 560 mm	XXX,-
JAY020011	JAY Balance, Fluid, Microclimatic Cover, Until Seat Width 500 mm	XXX,-
JAY020013	JAY Balance, Fluid, Microclimatic Cover, From Seat Width 560 mm	XXX,-
JAY020015	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm	XXX,-
JAY020017	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, From Seat Width 560 mm	XXX,-
JAY020018	JAY Xtreme Active	XXX,-
JAY020023	JAY Leg Support Kit For The Xtreme Active (2)	XXX,-
JAY020019	JAY J2 Cushion, Until Seat Width 460 mm Or Seat Depth 480 mm	XXX,-
JAY020021	JAY J2 Cushion, From Seat Width 480 mm Or Seat Depth 500 mm	XXX,-
JAY020020	JAY J2 Cushion, Deep Contour, Until Seat Width 480 mm	XXX,-
JAY020022	JAY J2 Cushion, Deep Contour, From Seat Width 500 mm	XXX,-

### Knee & Thigh Supports

€

QOP020342	Knee Support, Lateral, 120 x 120 mm, Swing Away, Right (3)	XXX,-
QOP020341	Knee Support, Lateral, 120 x 120 mm, Swing Away, Left (3)	XXX,-
QOP020344	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Right (4)	XXX,-
QOP020343	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Left (4)	XXX,-
QOP020347	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Right (4)	XXX,-
QOP020346	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Left (4)	XXX,-

(1) Only With Seat Width 560 mm

(2) Only For Cushions With Seat Depth 420 - 500 mm

(3) Not With QOP050015-17 SEDEO PRO, Legrest, Swing Away, QOP050080 SEDEO PRO, Hanger

Receivers Only, QOP090122 Transfer Bar

(4) Not With Seat To Floor Height 415 mm, Not With QOP040035 Side Guard, Wedge Shaped, Padded



## Powered Seating

### Armrest Types

€

QOP040019	Armrest, Single Post, Flip Back (1)	–
QOP040021	Armrest, Single Post, Removable (1)	XXX,-
QOP040123	Armrest, Reclining, Flip Back	XXX,-

### Arm Pads

€

QOP040023	Arm Pad, Integral Foam, 300 mm, Right	XXX,-
QOP040024	Arm Pad, Integral Foam, 300 mm, Left	XXX,-
QOP040025	Arm Pad, Integral Foam, 400 mm, Right	–
QOP040026	Arm Pad, Integral Foam, 400 mm, Left	–
QOP040027	Arm Pad, Integral Foam, 450 mm, Right (2)	XXX,-
QOP040028	Arm Pad, Integral Foam, 450 mm, Left (2)	XXX,-
QOP040029	Omit Arm Pad, Right	XXX,-
QOP040030	Omit Arm Pad, Left	XXX,-

### Comfort Arm Pads

QOP040204	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right	XXX,-
QOP040203	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Left	XXX,-

### Hemiplegia Arm Pads

QOP040031	Arm Pad, Hemiplegia, Integral Foam, Right (3)	XXX,-
QOP040032	Arm Pad, Hemiplegia, Integral Foam, Left (4)	XXX,-

### Side Guard Pads For Armrest

€

#### Side Guard Pads For Armrest, Single Post, Flip Back / Removable

QOP040033	Side Guard (5)	XXX,-
QOP040034	Side Guard, Integral Foam (6)	XXX,-
QOP040035	Side Guard, Wedge Shaped, Padded (7)	XXX,-

#### Side Guard Pads For Armrest, Reclining, Flip Back

QOP040036	Side Guard, Mounted On Seat Rail Bracket (8)	XXX,-
-----------	--	-------

### Centre Mount Legrest

€

QOP050086	Centre Mount Legrest, Tool Adjustable, With Flip Up And Angle-Adjustable Footboard (9) (10)	XXX,-
-----------	---	-------

#### Pad Options

QOP050087	Centre Calf Support, Integral Foam	XXX,-
QOP050092	Calf Pad, Integral Foam, Pair	XXX,-

### Power Centre Mount Legrest

€

QOP050076	Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension (10) (11)	XXX,-
QOP050071	Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension, Footplate To Floor (10) (11)	XXX,-

#### Footplates

QOP050072	Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	XXX,-
QOP050073	Footplate, Dual, Height- And Angle-Adjustable, Large, 310 x 170 mm	XXX,-

#### Pad Options

QOP050075	Calf Pad, Padded, Pair	XXX,-
-----------	------------------------	-------

(1) Not With Seat To Floor Height 415 mm

(2) With QOP040019 Flip Back Armrest, Seat Width +20 mm

(3) Not With QOP110052 Joystick Mount, Right, QOP090040 Tray Table, Swing Away, QOP090046

Tray Table, Half-Side, Swing Away, QOP090110 Rear View Mirror, Right

(4) Not With QOP110053 Joystick Mount, Left, QOP090040 Tray Table, Swing Away, QOP090046

Tray Table, Half-Side, Swing Away, QOP090118 Rear View Mirror, Left

(5) Not With QOP040123 Armrest, Reclining, Flip Back

(6) Seat Width -20 mm Not With QOP040123 Armrest, Reclining, Flip Back

(7) Seat Width -80 mm, Only With QOP040021 Armrest, Single Post, Removable

(8) Only With QOP040123 Armrest, Reclining, Flip Back

(9) In Combination With Lift Or Tilt, Choose QOP050087 Centre Calf Support, Integral Foam Or

QOP050092 Calf Pad, Integral Foam, Pair

(10) Not With RWD

(11) Not With QOP110036 Remote Joystick, VR2; QOP110037 Remote Joystick, VR2, 2 Actuators &

QOP110040 Remote Joystick, VR2, 2 Actuators And Lights. In Case Of R-Net: QOP110032 Actuator

Control & Lights Module, R-Net, ISM-6L

## Powered Seating

### SEDEO PRO Swing Away Legrest

€

QOP050015	SEDEO PRO, Legrest, 70°, Swing Away	XXX,-
QOP050016	SEDEO PRO, Legrest, Manual Elevating, Swing Away	XXX,-
QOP050017	SEDEO PRO, Legrest, Power Elevating, Swing Away	XXX,-
QOP050080	SEDEO PRO, Hanger Receivers Only	XXX,-
QOP050000	Omit Legrest	XXX,-

#### Footplates

QOP050088	Footplate, For Seat Width 385 - 400 mm	XXX,-
QOP050089	Footplate, For Seat Width 420 - 460 mm	XXX,-
QOP050090	Footplate, For Seat Width 480 - 520 mm	XXX,-
QOP050091	Footplate, For Seat Width 540 - 560 mm	XXX,-
QOP050081	Footplates, Coupled (1) (2)	XXX,-

#### Support Options

QOP050082	Calf Support, Integral Foam, Depth-Adjustable, Pair	XXX,-
QOP050083	Calf Support, Padded, Depth-Adjustable, Pair	XXX,-
QOP050084	Heel Loop	XXX,-
QOP050085	Calf Strap (2)	-
QOP050200	Knee Pad, Integral Foam	XXX,-

#### Amputee Support (3)

QOP050013	Amputee Support, Padded, Right	XXX,-
QOP050014	Amputee Support, Padded, Left	XXX,-

### Positioning Supports & Accessories

€

QOP090128	Positioning Belt, Carseat Buckle, Short, 50 mm x 1700 mm	-
QOP090129	Positioning Belt, Carseat Buckle, Long, 50 mm x 2700 mm	XXX,-
QOP090110	Rear View Mirror, Right (4)	XXX,-
QOP090118	Rear View Mirror, Left (4)	XXX,-
QOP090009	Crutch Holder (5)	XXX,-
QOP090119	Brackets For Bag On Backrest Frame (5)	XXX,-
QOP090120	Bag For Armrest (6)	XXX,-
QOP090122	Transfer Bar (7)	XXX,-
QOP030605	Abduction Pommel (8)	XXX,-

#### Tray Table (Mounted Opposite Joystick)

QOP090040	Tray Table, Width-Adjustable, Swing Away	XXX,-
QOP090046	Tray Table, Half-Side, Swing Away (9)	XXX,-

#### Tray Table Accessories

QOP090125	Partial Padding For Tray Table	XXX,-
QOP090126	Raised Edges, Acrylic For Tray Table	XXX,-
QOP090127	Elbow Pads, Padded, Pair For Tray Table	XXX,-

(1) Choose Also A Footplate

(2) Only For QOP050015 SEDEO PRO, Legrest, 70°, Swing Away

(3) Extra On A Set Of Swing Away Legrests, Attached Lose

(4) Not In Combination With Same Entry Side QOP090040 Tray Table, Swing Away, QOP090046 Tray Table, Half-Side, Swing Away

(5) Not With QOP030086 Powered Recline, 90° - 130°, QOP030087 Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear

(6) Not With QOP040123 Armrest, Reclining, Flip Back

(7) Not With QOP020342/QOP020341 Knee Support, Lateral, Swing Away, Right/Left &amp;

QOP050015-17 SEDEO PRO, Legrest, Swing Away &amp; QOP050080 SEDEO PRO, Hanger Receivers Only

(8) Not With (Power) Centre Mount Legrests

(9) Not With QOP090125-127 Tray Table Accessories

**Joystick / Electronics****Remote Joystick – VR2**

€

QOP110036	Remote Joystick, VR2	–
QOP110037	Remote Joystick, VR2, 2 Actuators	XXX,-
QOP110040	Remote Joystick, VR2, 2 Actuators And Lights	XXX,-

**Joystick Mounting**

€

QOP110050	Joystick Mounting, Fixed	–
QOP110056	Joystick Mounting, Parallel, Swing Away, Capsulated	XXX,-

**Joystick Mounting Location**

€

QOP110052	Joystick Mounting Location, Right	–
QOP110053	Joystick Mounting Location, Left	–

**Joystick Knobs (Optional)**

€

QOP090505	Joystick Knob, Golf Ball, Ø 42 mm	XXX,-
QOP090506	Joystick Knob, T-Bar, 82 mm Wide	XXX,-
QOP090507	Joystick Knob, Soft Ball, Ø 70 mm	XXX,-
QOP090508	Joystick Knob, Ball, Ø 35 mm	XXX,-
QOP090045	Joystick Knob, Allround, Ball, Ø 48 mm	XXX,-
QOP090509	Joystick Knob, Chin Cup, Ø 37 mm	XXX,-
QOP090512	Joystick Knob, Stick, 80 x 15 mm	XXX,-
QOP090511	Joystick Knob, Mushroom, Ø 45 mm	XXX,-

**Electronic Programmers – VR2**

€

QOP090065	Mobility Programmer With Dongle, VR2	XXX,-
-----------	--------------------------------------	-------

**Electronic Accessories Miscellaneous**

QOP090063	DTT Hand Held Programmer For VR2 With Connector Cable	XXX,-
-----------	---	-------

**Power Module – VR2**

€

QOP110002	Power Module, VR2, 90 A	XXX,-
-----------	-------------------------	-------

**Lighting Module – VR2**

€

QOP110030	Lighting Module, VR2 (1)	XXX,-
-----------	--------------------------	-------

**Attendant Controls**

€

QOP110021	Attendant Control, Dual	XXX,-
QOP110020	Joystick Mounted As Attendant Control (2)	XXX,-

**Attendant Mounting Location**

QOP110025	Attendant Mounting Location, Right	–
QOP110026	Attendant Mounting Location, Left	–

**4-Actuator Control Box (Direct Actuator Control)**

€

QOP110054	4-Actuator Control Box Attached To Armrest, Right	XXX,-
QOP110055	4-Actuator Control Box Attached To Armrest, Left	XXX,-

(1) If LED Light And Indicator Is Chosen, Then This Is Required

(2) Choose Also Remote Joystick

**Joystick / Electronics****Remote Joystick – R-Net** 

€

QOP110015	Remote Joystick, R-Net, LED, Actuators	–
QOP110016	Remote Joystick, R-Net, LED, Actuators And Lights	–
QOP110017	Remote Joystick, R-Net, Colour LCD	–
QOP110018	Remote Joystick, R-Net, Colour LCD, Heavy Duty	XXX,-
QOP110013	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	XXX,-
QOP110023	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	XXX,-
QOP110019	Omit Remote Joystick, R-Net	XXX,-

**Joystick Mounting** 

€

QOP110050	Joystick Mounting, Fixed	–
QOP110056	Joystick Mounting, Parallel, Swing Away, Capsulated	XXX,-

**Joystick Mounting Location**

€

QOP110052	Joystick Mounting Location, Right	–
QOP110053	Joystick Mounting Location, Left	–

**Display Language Selection**

€

QOP119000	German	QOP119004	Portuguese	QOP119009	Swedish	QOP119013	Slovak
QOP119001	English	QOP119005	Italian	QOP119010	Greek	QOP119014	Polish
QOP119002	French	QOP119006	Dutch	QOP119011	Finnish		
QOP119003	Spanish	QOP119008	Norwegian	QOP119012	Czech		

**Joystick Knobs** (Optional)

€

QOP090505	Joystick Knob, Golf Ball, Ø 42 mm	XXX,-
QOP090506	Joystick Knob, T-Bar, 82 mm Wide	XXX,-
QOP090507	Joystick Knob, Soft Ball, Ø 70 mm	XXX,-
QOP090508	Joystick Knob, Ball, Ø 35 mm	XXX,-
QOP090045	Joystick Knob, Allround, Ball Ø 48 mm	XXX,-
QOP090509	Joystick Knob, Chin Cup, Ø 37 mm	XXX,-
QOP090512	Joystick Knob, Stick, 80 x 15 mm	XXX,-
QOP090511	Joystick Knob, Mushroom, Ø 45 mm	XXX,-

**Electronic Programmers – R-Net** 

€

QOP090066	PC Programmer With Dongle, R-Net	XXX,-
-----------	----------------------------------	-------

**Electronic Accessories Miscellaneous**

QOP090064	DTT Hand Held Programmer For R-Net With Connector Cable	XXX,-
-----------	---	-------

**Power Modules – R-Net** 

€

QOP110003	Power Module, R-Net, 90 A (1)	XXX,-
QOP110004	Power Module, R-Net, 120 A	XXX,-

**Seating / Lighting Module Control – R-Net** 

€

QOP110032	Actuator Control & Lights Module, R-Net, ISM-6L	XXX,-
-----------	---	-------

**Attendant Controls** 

€

QOP110021	Attendant Control, Dual	XXX,-
QOP110020	Joystick Mounted As Attendant Control (1)	XXX,-

**Attendant Mounting Location**

QOP110025	Attendant Mounting Location, Right	–
QOP110026	Attendant Mounting Location, Left	–

(1) Choose Also Remote Joystick

## Joysticks / Electronics

### Table Top Controls

€

QOP110143	Table Top Control, LCD, Centre Tray Mounted	XXX,-
-----------	---	-------

### Tray Table Mounting Location

QOP110080	Tray Table Mounting Location, Left	–
QOP110081	Tray Table Mounting Location, Right	–

### #1 Ctrl+5 Switch Box Options

€

QOP110043	#1 Ctrl+5 Switch Box With Free Assignable Functions	XXX,-
-----------	---	-------

### #1 Ctrl+5 Switch Box Mounting Location

QOP110180	#1 Ctrl+5 Switch Box Mounting Location, Behind Joystick (1)	–
QOP110181	#1 Ctrl+5 Switch Box Mounting Location, Opposite Joystick	–

### #1 Ctrl+5 Switch Type

QOP110175	#1 Ctrl+5 Switch Type With Push Buttons	–
QOP110176	#1 Ctrl+5 Switch Type With Toggle Switches	–

### Bluetooth Mouse Modules / Electronics – R-Net

€

QOP090068	Bluetooth Mouse Module, R-Net	XXX,-
QOP090141	2nd Bluetooth Mouse Module, R-Net (Android And Windows)	XXX,-
QOP090142	iOS Bluetooth Mouse Module, R-Net (Apple, iOS)	XXX,-

## Special Controls

### OMNI Input Control Module – R-Net <sup>(2)</sup>

€

QOP110101	OMNI Input Control Module, R-Net (3)	XXX,-
QOP110102	OMNI Input Control Module, R-Net, Scanning Mode	XXX,-
QOP110103	OMNI Input Control Module, R-Net, Scanning Mode And Traction Control System (GYRO Technology)	XXX,-
QOP110008	Programming Service Not Discountable	XXX,-

### OMNI Module Mount <sup>(4)</sup>

€

QOP110126	OMNI Module, Fixed Mount, Right	–
QOP110127	OMNI Module, Fixed Mount, Left	–
QOP110128	OMNI Module, Swing Away Mount, Right	XXX,-
QOP110129	OMNI Module, Swing Away Mount, Left	XXX,-

(1) Not Possible With QOP110056 Joystick Mounting, Parallel Swing Away, Capsulated (Here We Can't Swing Away)

(2) Choose OMNI Mount And Input Device in Sections Below And Choose Switches in Section E

(3) Display With Built In Infrared & Armrest Mount, Required For Specialty Input Devices

(4) If Used With Standard Joystick Selection, It Must Be Mounted Opposite Side

**Hand Control (1)**

€

**Proportional Input Devices**

QOP110301	Hand Control, Switch-It MicroGuide	XXX,-
QOP110302	Hand Control, Switch-It MicroPilot	XXX,-
QOP110303	Hand Control, Switch-It TouchDrive 2	XXX,-
QOP110304	Hand Control, Switch-It VersaGuide	XXX,-
QOP110305	Hand Control, Switch-It VersaGuide EZ (2)	XXX,-
QOP110306	Hand Control, Allround Joystick	XXX,-
QOP110307	Hand Control, Allround Joystick Light (3)	XXX,-
QOP110308	Hand Control, Allround Joystick Heavy Duty	XXX,-
QOP110309	Hand Control, Allround Joystick, R-Net Bus Enabled (4)	XXX,-
QOP110310	Hand Control, Allround Joystick Light, R-Net Bus Enabled (4)	XXX,-
QOP110311	Hand Control, Allround Joystick Heavy Duty, R-Net Bus Enabled (4)	XXX,-

**Proportional Input Devices, Price Inclusive OMNI Module And Traction Control System (GYRO Technology)**

QOP110044	Hand Control, Switch-It MicroGuide	XXX,-
QOP110045	Hand Control, Switch-It MicroPilot	XXX,-
QOP110046	Hand Control, Switch-It TouchDrive 2	XXX,-
QOP110047	Hand Control, Switch-It VersaGuide	XXX,-
QOP110048	Hand Control, Switch-It VersaGuide EZ (2)	XXX,-
QOP110049	Hand Control, Allround Joystick	XXX,-
QOP110059	Hand Control, Allround Joystick Light (3)	XXX,-
QOP110060	Hand Control, Allround Joystick Heavy Duty	XXX,-

**Hand Control Mounting Platforms**

€

QOP110342	Hand Control Mounting, Link-It To The Space Disk, Hand Tray (5)	XXX,-
-----------	---	-------

**Hand Control Mounting **

QOP110350	Hand Control Mounting, Link-It To The Armrest, Fixed	XXX,-
QOP110353	Hand Control Mounting, Link-It To The Armrest, Swing Away	XXX,-

**Hand Control Mounting Location**

QOP110362	Hand Control Mounting Location, Right	XXX,-
QOP110361	Hand Control Mounting Location, Left	XXX,-

**Head Control (6) **

€

**Proportional Input Devices**

QOP110370	Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	XXX,-
QOP110371	Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads	XXX,-
QOP110372	Head Control, Switch-It RIM, 3-Way Proportional	XXX,-

**Proportional Input Devices, Price Inclusive OMNI Module And Traction Control System (GYRO Technology)**

QOP110061	Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	XXX,-
QOP110062	Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads	XXX,-
QOP110063	Head Control, Switch-It RIM, 3-Way Proportional	XXX,-

(1) For Switch Selection See Section E

(2) Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide

(3) Only Use With Default Joystick Knob Or Joystick Knob, Allround, Ball Ø 48 mm

(4) OMNI Module Not Required

(5) Only For Micro Series Joysticks

(6) For Switch Selection See Section E. Choose A Whitmyer Headrest Mount On Previous Pages

**Switched Input Devices**

QOP110401	Head Control, 3-Switch Head Array, Wings, Short, Curved	XXX,-
QOP110402	Head Control, 3-Switch Head Array, Wings, Short, Straight	XXX,-
QOP110403	Head Control, 3-Switch Head Array, Wings, Long, Curved	XXX,-
QOP110404	Head Control, 3-Switch Head Array, Wings, Long, Straight	XXX,-
QOP110405	Head Control, 4-Switch Head Array, Wings, Long, Straight	XXX,-
QOP110406	Head Control, 5-Switch Head Array, Wings, Long, Straight	XXX,-

**Switched Input Devices, Price Inclusive OMNI Module And Traction Control System (GYRO Technology)**

QOP110064	Head Control, 3-Switch Head Array, Wings, Short, Curved	XXX,-
QOP110065	Head Control, 3-Switch Head Array, Wings, Short, Straight	XXX,-
QOP110066	Head Control, 3-Switch Head Array, Wings, Long, Curved	XXX,-
QOP110067	Head Control, 3-Switch Head Array, Wings, Long, Straight	XXX,-
QOP110068	Head Control, 4-Switch Head Array, Wings, Long, Straight	XXX,-
QOP110069	Head Control, 5-Switch Head Array, Wings, Long, Straight	XXX,-

**Chin Control** 

€

**Proportional Input Devices**

QOP110470	Chin Control, Switch-It MicroGuide	XXX,-
QOP110471	Chin Control, Switch-It MicroPilot	XXX,-
QOP110472	Chin Control, Switch-It VersaGuide	XXX,-
QOP110473	Chin Control, Switch-It VersaGuide EZ	XXX,-
QOP110474	Chin Control, Allround Joystick	XXX,-
QOP110475	Chin Control, Allround Joystick Light	XXX,-
QOP110476	Chin Control, Allround Joystick, R-Net Bus Enabled (1)	XXX,-
QOP110477	Chin Control, Allround Joystick Light, R-Net Bus Enabled (1)	XXX,-

**Proportional Input Devices, Price Inclusive OMNI Module And Traction Control System (GYRO Technology)**

QOP110070	Chin Control, Switch-It MicroGuide	XXX,-
QOP110071	Chin Control, Switch-It MicroPilot	XXX,-
QOP110072	Chin Control, Switch-It VersaGuide	XXX,-
QOP110073	Chin Control, Switch-It VersaGuide EZ	XXX,-
QOP110074	Chin Control, Allround Joystick	XXX,-
QOP110075	Chin Control, Allround Joystick Light	XXX,-

**Chin Control Mounting** 

QOP110480	Chin Control Mounting, Link-It To Chin Harness, Small (2)	XXX,-
QOP110481	Chin Control Mounting, Link-It To Chin Harness, Large (2)	XXX,-
QOP110483	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Right	XXX,-
QOP110484	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Left	XXX,-
QOP110485	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Right	XXX,-
QOP110486	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Left	XXX,-
QOP110487	Chin Control Mounting, Link-It To Headrest, Swing Away, Left	XXX,-
QOP110488	Chin Control Mounting, Link-It To Headrest, Swing Away, Right	XXX,-

**Link-It Extension Kit**

QOP110590	Link-It Extension Kit, Base Module, 100 mm Extension Tube, 5 Cable Clips	XXX,-
QOP110591	Link-It Extension Kit, Base Module, 150 mm Extension Tube, 5 Cable Clips	XXX,-

**Sip 'N Puff** (3)

€

QOP110510	Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest (4)	XXX,-
-----------	---	-------

**Sip 'N Puff, Price Inclusive OMNI Module And Traction Control System (GYRO Technology)**

QOP110076	Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest (4)	XXX,-
-----------	---	-------

(1) OMNI Module Not Required

(2) If Switches Required Only Allround Joystick With Satellite Switches Available

(3) For Switch Selection See Section E

(4) Headrest To Be Selected

## Special Controls

E – Switches (1) 

€

## #1 Switch Selection

QOP110520	#1 Buddy Button – Black	XXX,-
QOP110521	#1 Buddy Button – Red	XXX,-
QOP110522	#1 Buddy Button – Green	XXX,-
QOP110523	#1 Egg Switch – Black	XXX,-
QOP110524	#1 MicroLite Switch – Black	XXX,-
QOP110526	#1 Piko Button – Black	XXX,-
QOP110527	#1 Piko Button – Blue	XXX,-
QOP110528	#1 Piko Button – Green	XXX,-
QOP110529	#1 Piko Button – Red	XXX,-
QOP110530	#1 Piko Button – Yellow	XXX,-
QOP110534	#1 Twister Switch – Black	XXX,-
QOP110535	#1 Twister Switch – Green	XXX,-
QOP110536	#1 Twister Switch – Red	XXX,-

#1 Switch Mounting 

QOP110560	#1 Switch Mounting, Link-It To The Hand Control (2)	XXX,-
QOP110566	#1 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (3)	XXX,-
QOP110561	#1 Switch Mounting, Link-It To The Headrest	XXX,-
QOP110563	#1 Switch Mounting, Link-It To The Seat Rail (4)	XXX,-
QOP110565	#1 Switch Mounting, Link-It To Chin Boom (5)	XXX,-
QOP110569	#1 Omit Switch Mounting (6)	XXX,-

## #1 Switch Mounting Location

QOP110572	#1 Switch Mounting Location, Right	XXX,-
QOP110571	#1 Switch Mounting Location, Left	XXX,-

## #1 Satellite Switches (7)

QOP110537	#1 Satellite Switch – Black	XXX,-
QOP110538	#1 Satellite Switch – Green	XXX,-
QOP110539	#1 Satellite Switch – Red	XXX,-

## #2 Switch Selection

QOP110083	#2 Buddy Button – Black	XXX,-
QOP110086	#2 Buddy Button – Red	XXX,-
QOP110087	#2 Buddy Button – Green	XXX,-
QOP110088	#2 Egg Switch – Black	XXX,-
QOP110089	#2 MicroLite Switch – Black	XXX,-
QOP110090	#2 Piko Button – Black	XXX,-
QOP110091	#2 Piko Button – Blue	XXX,-
QOP110092	#2 Piko Button – Green	XXX,-
QOP110093	#2 Piko Button – Red	XXX,-
QOP110094	#2 Piko Button – Yellow	XXX,-
QOP110095	#2 Twister Switch – Black	XXX,-
QOP110096	#2 Twister Switch – Green	XXX,-
QOP110097	#2 Twister Switch – Red	XXX,-

## #2 Switch Mounting

QOP110098	#2 Switch Mounting, Link-It To The Hand Control (1)	XXX,-
QOP110099	#2 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (2)	XXX,-
QOP110100	#2 Switch Mounting, Link-It To The Headrest	XXX,-
QOP110104	#2 Switch Mounting, Link-It To The Seat Rail (3)	XXX,-
QOP110105	#2 Switch Mounting, Link-It To Chin Boom (4)	XXX,-
QOP110106	#2 Omit Switch Mounting (5)	XXX,-

(1) Select 1 Or 2 Switches And Mounting Positions Accordingly. Single Function For Special Input Devices

(2) Not With HD Allround Joystick Or Space Disk Mount

(3) Armrest Rail Mount

(4) Knee, If 2 Switches Selected Only Left And Right Mount

(5) No Left / Right Switch Selection Required

(6) Switches Attached Lose

(7) Only For Allround Joystick, Maximum 2 Selectable

(8) Used In Combination With Assignable Button Programming

(9) Attached Lose, When Ordered. No Link-It Fixation Possible



## Special Controls

### Multi-Function Switches (8)

QOP110550	Mini 2 Button (Stereo Plug) (9)	XXX,-
QOP110551	Mini 2 Button (Dual Mono Plug) (9)	XXX,-
QOP110552	Mini 4 Button InLine (DB9 Connector) (9)	XXX,-
QOP110554	4-Way Toggle Switched Joystick	XXX,-

### Accessories

QOP110580	Mono Extension Cable (3.5 mm Jack Socket) (2 m)	XXX,-
QOP110581	Stereo Extension Cable (3.5 mm Jack Socket) (2 m)	XXX,-
QOP110582	Extension Cable DB9 (1 m)	XXX,-
QOP110583	Y-Cable, 2 Female Mono To 1 Male Stereo	XXX,-

## Safety Devices

€

QOP090147	Remote Stop Switch (30 m Activation Range, Fit Via OBC Port) (10)	XXX,-
-----------	---	-------

**Drive Wheels & Castor Wheels & Mudguards & Kerb Climber**

◀ BACK



// QOP070069



// QOP080213



// QOP010007



// QOP010008



// QOP90038

**LED Lights & Chargers**

◀ BACK



// QOP110209



// QOP110209



// QOP160101



// QOP160104

**SEDEO PRO Back Cushions & Interfaces**

◀ BACK



// QOP030109



// QOP030113



// QOP030114



// QOP030115



// QOP030116



// QOP030117

## Headrest Pads & Headrest Mounts

◀ BACK




// QOP030430 & QOP030431



// WHT030401-WHT030403



// WHT030404-WHT030406



// QOP030433



// WHT030408



// WHT030409

## SEDEO Seat Cushions

◀ BACK



// QOP020073



// QOP020075



// QOP020076



// QOP020077



// QOP020078

**JAY Seat Cushions**

◀ **BACK**



// JAY020005



// JAY020003



// JAY020006



// JAY020010



// JAY020011



// JAY020018




// JAY020019



// JAY020020

**Armrest Types & Arm Pads**

◀ **BACK**



// QOP040019



// QOP040021



// QOP040123



// QPP040204



// QOP040027



// QOP040023

**Power Centre Mount Legrest & Centre Mount Legrest**

◀ **BACK**




// QOP050076 & QOP050071



// QOP050086

**SEDEO PRO Swing Away Legrest & Support Options**

◀ BACK



// QOP050015



// QOP050016



// QOP050082

**Remote Joystick– R-Net & Joystick Mounting**

◀ BACK



// QOP110017



// QOP110013



// QOP110018



// QOP110056

**Remote Joystick – VR2**

◀ BACK



// QOP110036



// QOP110037

**Joystick Knobs**

◀ BACK



// QOP090505



// QOP090045



// QOP090507



// QOP090508



// QOP090509



// QOP090512



// QOP090511



// QOP090506

**Electronic Programmers – VR2**

◀ BACK



// QOP090063

**Electronic Programmers – R-Net**

◀ BACK



// QOP090066



// QOP090064

**Power Modules – R-Net & Seating / Lighting Module Control – R-Net**

◀ BACK



// QOP110004

**Attendant Controls**

◀ BACK



// QOP110021

**#1 Ctrl+5- Switch Box Options & OMNI Input Control Module – R-Net**

◀ BACK



// QOP110175



// QOP110176



// QOP090068



// QOP110101

## Head Control

[← BACK](#)


// QOP110401 - QOP110406



// QOP110370 // QOP110371

## Chin Control & Hand Control

[← BACK](#)


// QOP110470 // QOP110301



// QOP110471 // QOP110302

// QOP110472 // QOP110473 //  
QOP110304 // QOP110305// QOP110306 - QOP110309  
// QOP110474 - QOP110477// QOP110475 // QOP110477  
// QOP110307 // QOP110310

## Table Top Controls

[← BACK](#)


// QOP110143



Switches

◀ BACK



// QOP110521



// QOP110523



// QOP110524



// QOP110526 - QOP110530



// QOP110537 - QPP110539



// QOP110534 - QOP110536



// QOP110551



// QOP110552

Switch Mounting

◀ BACK



// QOP110526 // QOP110566



// QOP110566 // QOP110526  
// QOP110529



// QOP110560 // QOP110572  
// QOP110013 // QOP110056  
// QOP110529 // QOP110526



// QOP110563



**Hand Control Mounting**

◀ **BACK**



// QOP110353 // QOP110302  
// QOP110362



// QOP110353 // QOP110302  
// QOP110362 // QOP110342



// QOP110302



// QOP110353 // QOP110361  
// QOP110524 // QOP110526  
// QOP110560



// QOP110353 // QOP110306  
// QOP110362 // QOP110538  
// QOP110539

**Chin Control Mounting**

◀ **BACK**



// QOP110480 // QOP110481  
// QOP110470 // QOP110471



// QOP110483 // QOP110484



// QOP110483 // QOP110484



// QOP110485 // QOP110486



// QOP110485 // QOP110474



// QOP110487 // QOP110471  
// QOP110526 // QOP110529  
// QOP110561 // QOP110572  
// QOP110571



// QOP110487 // QOP110474

## Comments For Customisations



## Q500 M/F/R with SEDEO PRO

€

<b>Base Price</b>	<b>00,000.00</b>
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
<b>Total Price</b>	

Date

Signature

Day      Month      Year


**SUNRISE**  
 MEDICAL

\*Sunrise Medical Limited

Thorns Road  
 Brierley Hill  
 DY5 2LD/West Midlands  
 Phone: + 44 (0) 845 605 6688  
 Fax: + 44 (0) 845 605 6689  
 E-Mail: Orders@sunmed.co.uk  
 www.SunriseMedical.eu

Sunrise Medical GmbH

Kahlbachring 2-4  
 D-69254 Malsch/HD  
 Phone: +49 (0) 72 53 / 9 80-0  
 Fax: +49 (0) 72 53 / 9 80-222  
 E-Mail: export@sunrisemedical.de  
 www.SunriseMedical.eu