

FORA TD-3124

BLOOD PRESSURE MONITORING SYSTEM



FEATURES

- Clever inflation technology
- 100 sets of memory (TD-3124B only)
- IHB detection
- One button and large LCD display design

SPECIFICATIONS

Accuracy of Pressure	±3mmHg or ±2% of reading
Accuracy of Heart Rate	±4% of reading
Measurement Range	Systolic: 50 ~ 250 mmHg / Diastolic: 30 ~ 180 mmHg Pulse rate: 40 ~ 199 / minute
Power Source	4 x AA (TD-3124A) / 4 x AA or 6V DC adaptor (optional for TD-3124B)
Dimension	130.5 (L) x 90 (W) x 55 (H) mm
Weight	210 g (without battery)
Operating Condition	10°C ~ 40°C, below 85% R.H.
Power Saving	180 seconds
Unit per Carton	30
Net Weight per Carton	10.8 kg
Gross Weight per Carton	11.8 kg
Cuff	Wide Range 24-43cm

