


SURFLO® Winged Infusion Sets

Our SURFLO Wing design minimizes the risk of vessel damage and enhances ease of penetration.

Compare our wings versus the competition and see why SURFLO is a highly respected brand.


Terumo's soft, flexible interlocking wings readily conform to the body's contour. The wing is angled for ease of positioning of the needle during venipuncture, for less risk of vessel damage or occlusion of the lumen.

The precisely honed needle with ultra-thin wall design assures ease of penetration and continued access with a single insertion.



These products contain no components made of natural rubber latex.

Surflo® Winged Infusion Sets



Product Code Description

3½" Tubing

Packed: 100/box, 5 boxes/case (500)

Product Code	Description
SV*19BLS	19G x ¾" TW
SV*21BLS	21G x ¾" UTW
SV*23BLS	23G x ¾" UTW
SV*25BLS	25G x ¾" TW

RW - Regular Wall
 TW - Thin Wall
 UTW - Ultra Thin Wall

Product Code Description

8" Tubing

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
SV*25EL	25G x ½" TW
SV*27EL	27G x ½" RW

Product Code Description

12" Tubing

Packed: 50/box, 10 boxes/case (500)

Product Code	Description
SV-18BLK	18G x ¾" TW
SV-19BLK	19G x ¾" TW
SV-21BLK	21G x ¾" UTW
SV-22BLK	22G x ¾" UTW
SV-23BLK	23G x ¾" UTW
SV-25BLK	25G x ¾" TW

Call your Terumo representative at:
1.800.888.3786
 or visit www.TerumoTMP.com
 and click on Find Your Local Terumo Rep.

TERUMO and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. SURFLO is a trademark owned by Terumo Corporation, Elkton, MD. PM-00689

