



E | R10

**PUSH SMARTER, NOT HARDER,
WITH EVERY STEP.**



**QUICK
INSTALLATION
GUIDE**



E | R10

The power assist device designed specifically to support the caregivers.



LIGHTWEIGHT

At just under 6.9 kgs. the R10 is one of the lightest power assists on the market.

PERFORMANCE

No compromise on performance. The R10 achieves a 12 km range, climbs up to 8° slopes, and handles a max user weight of 130 kg.

RELIABILITY

A solid, long lasting tire, lithium-ion battery, and a 3 year motor warranty promising quality and reliability.



TURN SLOPES FLAT AND HEAVY LOADS WEIGHTLESS.

Inclines can be especially challenging for wheelchair attendants. But with the R10's **Constant Speed Assistance** mode, you can get a continuous speed at the press of a button – even when tackling inclines of up to 8 degrees!



THERE'S NO NEED TO PRESS A BUTTON.

Upon detecting motion, **Automatic Assistance** provides just the right amount of propulsion to match your pace, reducing the effort needed to push. The faster you walk, the more assistance you receive. Plus, because you're only supplied power when you need it, your battery lasts longer too.



DISTANCE DETECTION SENSOR.

The R10's distance detection sensor automatically pauses the device if it doesn't detect an attendant gripping the push handle. This ensures the wheelchair only moves when someone is in control, adding an extra layer of safety for both you and the wheelchair user.

Compatible with almost all common manual wheelchairs with the following dimensions:

Box frame or similar frame type with round tubes diameter: 20 - 28 mm.

Seat width: 380 - 520 mm

Rear seat height: 400 - 580 mm

Height between the box frame and the ground: 180 - 220 mm



Quick Release of Control Unit

Ensures easy and safe transport without disconnecting cables.



Modern Back Pack

Compact and easy transport adapts to nearly any wheelchair.



Simple Attachment and Removal

Attaches easily for docking via small pin. Low added weight & simple design.



Easy & Convenient Folding

The R10 can remain in place when folded.



Simple Control

Thumb throttle and three push-button speed settings.

R10 Specifications

Weight Capacity	130 kg
Range*	up to 12 km
Max Slope	8°
Weight (including battery)	6.9 kg
Weight (battery pack)	1.2 kg
Adjustable Max Speed	1.5 kmh 3.0 kmh 4.5 kmh
Motor	250 Watt 36V / 15 Amp Brushless Motor
Battery	Lithium Ion 5 Ah / 36 V
Charging Time	Approx. 2.5 - 3 hours
Wheel Size	6"
Modes	Automatic Assistance Constant Speed (Cruise Control)
Warranty	3 year Motor 2 year Battery

See order form as specifications may vary based upon configuration.

* Depending on user weight and topography



Sunrise Medical Australia Pty Ltd
11 Daniel Street,
Wetherill Park, NSW 2164
Phone: (02) 9678 6600
Email: enquiries@sunrisemedical.com.au
Website: www.sunrisemedical.com.au

Manufacturer: Guangzhou ZUBU Medical Technology Co. Ltd, 11th Floor, Building 5, No. 2 Fangda Road, Huangpu Distract, Guangzhou/China. For further information on the full specification, options, and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

EMPOWER YOUR LIFE

