

TANITA MC780MA

- GMON
 - Create account
 - Save History
 - Show the result in different format
- eg, MC980MA format
- Data can output as excel format (which is user-friendly for research use)
- for data analyzer

A	B	C	D	E	F	G	H	I	J	K
CUSTOMERID	SEX	PATNR	DATETIME	MODEL	FLAG	CONDITION	AGE	BTYP	HEIGHT	WEIGHT
1	000	M	01/06/2018 15:50	MC-980		2	34.00	0	160.00	69.10
2	0000	F	09/10/2018 11:00	MC-780		2	24.00	0	159.00	52.80
3	0012345	F	16 21/01/2019 16:41	MC-780		2	26.00	0	160.00	67.40
4	0012345	F	16 24/10/2019 10:37	MC-780		2	27.00	0	160.00	66.70
5	0012345	F	16 24/10/2019 10:45	MC-780		2	27.00	0	160.00	66.80
6	0012345	F	16 15/03/2019 12:58	MC-780		2	26.00	0	160.00	67.60
7	0012345	F	6 21/05/2018 10:40	MC-780		2	25.00	0	160.00	68.50
8	1	F	6 05/09/2018 12:58	MC-980		2	26.00	0	160.00	68.10
9	1	F	6 02/10/2018 12:55	MC-980		2	26.00	0	160.00	67.60
10	1	F	6 02/10/2018 12:56	MC-980		2	26.00	0	160.00	67.60
11	1	F	6 02/10/2018 12:58	MC-980		2	26.00	0	160.00	67.60
12	1	F	6 02/10/2018 14:15	MC-980		2	26.00	0	160.00	68.00
13	1	F	6 02/10/2018 14:16	MC-980		2	26.00	0	160.00	68.00
14	1	F	6 30/10/2018 12:52	MC-980		2	26.00	0	160.00	68.00
15	1	F	6 03/01/2019 12:46	MC-980		2	26.00	0	160.00	67.10
16	1	F	6 22/02/2019 13:05	MC-980		2	26.00	0	160.00	66.00
17	1	F	6 04/03/2019 12:13	MC-980		2	26.00	0	160.00	66.80
18	1	F	6 04/03/2019 12:17	MC-980		2	26.00	0	160.00	66.80
19	1	F	6 04/03/2019 12:18	MC-980		2	26.00	0	160.00	66.80
20	1	F	6 29/03/2019 11:15	MC-980		2	26.00	0	160.00	65.80
21	1	F	6 29/03/2019 11:16	MC-980		2	26.00	0	160.00	65.80
22	1	F	6 29/03/2019 11:18	MC-980		2	26.00	0	160.00	65.80
23	1	F	6 29/03/2019 13:03	MC-980		2	26.00	0	160.00	65.80
24	1	F	6 23/04/2019 12:10	MC-980		2	26.00	0	160.00	65.30
25	1	F	6 29/04/2019 10:29	MC-980		2	26.00	0	160.00	65.80
26	1	F	6 29/04/2019 10:30	MC-980		2	26.00	0	160.00	65.80
27	1	F	6 29/04/2019 10:30	MC-980		2	26.00	0	160.00	65.80

資料庫 工具 說明

TANITA PRO

搜尋: 在 用戶編號-ID 群組: 所有群組

英

用戶編號-ID	姓	名	出生日期	SEX	城市
000			29/4/1984	M	
0000	Chan	Jessuca	16/9/1994	F	
0001		ciel	29/5/1997	F	
0012345		wing	9/6/1992	F	
1	w		9/6/1992	F	
123		STEP	22/10/1996	F	
1233			24/9/1988	M	
12334455			3/6/1973	F	
1234		elaine	8/7/1960	F	
12345	Dr Liu		25/8/1966	M	
123456	cc	hhh	7/6/1992	F	
1234567	chan	shuk man	30/7/1989	F	
12345678	so		27/10/1969	M	
1234678		cheung	18/7/1974	F	
1234777			17/1/1981	M	
123www	Peter		1/1/1995	M	
2182019		Mo	4/8/1986	M	

軟體資訊

評估風險記錄: 1
週期: 1/6/2018 到 1/6/2018 (0 個月)

© Medizin & Service GmbH - 版本 3.4.2

健康監控系統 - 版本 3.4.2 - Tanita Demo

資料 工具 說明

wing

身體測量 身體數值 進生數值 TANITA 家庭型體組成計 TANITA 專業型體組成計 活動量計 理想目標值 用戶資料

ID: 0012345
wing
年齡: 27.5歲

測量結果 分析

體組成分析 MC-980

測量狀態: 一般

體型: 一般

體重 [kg]: 64.0

身高 [cm]: 160

體脂率 [%]: 12.3

細胞外液 (ECW) [kg]: 18.2

BMI=25.0
細胞外液比 (ECW/TBW)=40.3%
四肢骨骼肌 (SMM)=7.73 kg/m²
骨骼肌量 (SMM)=24.4 kg (38.1%)

部位分析

	右腿	左腿	右
體脂率 [%]:	34.5	34.8	32.6
脂肪量 [kg]:	4.50	4.50	1.00
脂肪量 (FFM) [kg]:	8.60	8.40	2.00
肌肉量 (MM) [kg]:	8.10	7.90	1.90
電阻值 [ohm]:	224	229	369

列印報告

選擇印表機: Kyocera FS-C8520MFP X0

選項

MC-980 格式

預防檢查

體組成分析

MC-100 格式

MC-780 格式

MC-980 格式

BC-420 格式

綜合報告

體組成報告

特殊評估

體重

複製: 1

列印

儲存檔案

預覽

E-Mail 傳送 PDF 檔

取消

狀態: 19/12/2019 11:23:30

> 新測量

> 儲存

> 健康風險

> 風險值表

> 報告