


MC-780



Date / Time	Age 30 Standard ☒ Athletic ☐
17/5/2022 15:37	
Customer-ID / Name	Height 171.0 cm Male ☒ Female ☐
sample11	Clothes Weight (PT) 0.2 kg

Details





MC-780	Result	Desirable	Target	Diff
 Body Weight	67.7 kg	54.1-73.1		
 Fat %	22.9 %	11.0-22.0		
Fat Mass	15.5 kg	6.5-14.7		
FFM	52.2 kg			
 Muscle Mass	49.5 kg	47.5-55.7		
 Bone Mass	2.7 kg			
 BMI	23.2	18.5-25.0		
SMM	26.3 kg	Detailed Explanation on www.mc780.gmon.eu		
 Metabolic Age	36			
				



☐ **TBW - ECW** (Extracellular Water) - **ICW** (Intracellular Water) - **BMR** - **VFR**



❑ Segmental Analysis

☐ **Balance**

	<i>BIA Information</i>			
	5kHz	50kHz	250kHz	Phase Angle
H-L	782.2 / 46.0	672.1 / 78.0	599.4 / 65.3	6.6°
RL	324.9 / 16.5	281.2 / 30.2	253.4 / 21.6	6.1°
LL	332.6 / 20.0	285.1 / 31.6	256.5 / 22.8	6.3°
RH	407.5 / 24.9	347.6 / 42.8	304.7 / 38.3	7.0°
LH	417.8 / 22.2	360.3 / 42.4	317.3 / 39.7	6.7°
L-L	656.3 / 36.3	564.2 / 63.4	505.7 / 44.9	6.4°
RH-RL	762.8 / 44.3	655.0 / 76.2	582.3 / 62.2	6.6°