

eTSE

Toilet stand-up aid

The eTSE makes it easier to stand up thanks to a gentle tilting movement and can be operated continuously using a remote control. In combination with our TSE 120 or TSE 150, it is the perfect addition to the classic toilet seat raiser.

- Frameless, compact and space-saving design – does not require any additional space to the existing toilet seat raiser
- Continuously adjustable lifting and simultaneous tilting for particularly gentle and safe standing up
- Battery-powered, no permanent power supply required
- Controlled by remote control on the armrest
- Our TSE 120 or TSE 150 is required for use (not included)
- Simple, tool-free installation



TECHNICAL SPECIFICATIONS

Article number	520300000
Dimensions (L x W x H)	8,7 cm x 40,6 cm x 15,8 cm
Vertical lift	up to 95 mm
Tilt	up to 20°
Height adjustment	stepless
Total weight	4,25 kg
Max. load	150 kg